



NATtrol™

IVD MOLECULAR CONTROLS

FOR THOSE WHO MUST GET IT RIGHT

ZeptoMetrix is your end-to-end partner for high-quality infectious disease controls for your IVD diagnostic assays. Our ZeptoMetrix NATtrol™ is the most trusted brand of third-party molecular diagnostic QC for infectious disease testing. Using NATtrol Molecular Controls on a consistent basis helps to monitor testing performance and provides additional confidence in the quality and reliability of clinical laboratory operations.

What really matters to your lab: Refrigerated ready-to-use + Non-infectious + Whole microorganisms + Mimics a true clinical specimen










- Can be used across various molecular testing platforms and assays
- Available in different formats (single or multiple analytes)
- Full process controls that monitor extraction and amplification steps
- Train and monitor laboratory personnel
- Evaluate lot-to-lot consistency of test kits and assay reagents







ZEPTOMETRIX.COM +1.800.274.5487

ZeptoMetrix[®]
an antylia scientific company








Respiratory

Part Number	Product Description	Use*	Platform(s)**
NATFLUAB-6C-IVD	Influenza A/B Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATCXVA9-6C-IVD	Influenza/RSV Negative Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATCV9-6C-IVD	Influenza/RSV/SARS-CoV-2 Negative Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®, Roche cobas® Liat®
NATFRC-6C-IVD	Influenza/RSV/SARS-CoV-2 Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®, Roche cobas® Liat®
NATRSV-6C-IVD	Respiratory Syncytial Virus Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATSARS(COV2)-NEG-IVD	SARS-CoV-2 Negative Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®, Real-time PCR
NATSARS(COV2)-NEG1-IVD	SARS-CoV-2 Negative Control (6 x 1 mL)	CE IVD 	Cepheid® GeneXpert®, Real-time PCR
NATSARS(COV2)-ERC-IVD	SARS-CoV-2 Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®, Real-time PCR
NATSARS(COV2)-ERC1-IVD	SARS-CoV-2 Positive Control (6 x 1 mL)	CE IVD 	Cepheid® GeneXpert®, Real-time PCR






Gastrointestinal (GI)

Part Number	Product Description	Use*	Platform(s)**
NATCSO-6MC-IVD	C. difficile Negative Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATCDI-6MC-IVD	C. difficile Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATNOV-6MC-IVD	Norovirus GI/GII Positive Control 6 X 0.5 mL	CE IVD 	Cepheid® GeneXpert®
NATROTA-6MC-IVD	Norovirus Negative Control 6 X 0.5 mL	CE IVD 	Cepheid® GeneXpert®

Women's Health & STI

Part Number	Product Description	Use*	Platform(s)**
NATCT(434)-6MC	C. trachomatis Positive Control (6 x 1 mL)	IVD	Cepheid® GeneXpert®, binx® io
NATCT(434)-6MC-IVD	C. trachomatis Positive Control Pack (6 x 1.25 mL)	CE IVD 	Cepheid® GeneXpert®, Hologic® Panther Fusion®
NATCTNGNEG-6MC	C. trachomatis/N. gonorrhoeae Negative Control (6 x 1 mL)	IVD	Cepheid® GeneXpert®
NATCT-NGNEG-6MC-IVD	C. trachomatis/N. gonorrhoeae Negative Control (6 x 1.25 mL)	CE IVD 	Cepheid® GeneXpert®, Hologic® Panther Fusion®
NATCTNGTV-NEG-IVD	CT/NG/TV Negative Control (6 x 1.0 mL)	IVD	Visby Medical™
NATCTNGTV-POS-IVD	CT/NG/TV Positive Control (6 x 1.0 mL)	IVD	Visby Medical™
NATLAC-6MC-IVD	GBS Negative Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATSAG-6MC-IVD	GBS Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATNG-6MC	N. gonorrhoeae Positive Control (6 x 1 mL)	IVD	Cepheid® GeneXpert®, binx® io
NATNG-6MC-IVD	N. gonorrhoeae Positive Control (6 x 1.25 mL)	CE IVD 	Cepheid® GeneXpert®
NATTVNEG-6MC-IVD	T. vaginalis Negative Control (6 x 1.2 mL)	CE IVD 	Cepheid® GeneXpert®
NATTVPOS-6MC-IVD	T. vaginalis Positive Control (6 x 1.2 mL)	CE IVD 	Cepheid® GeneXpert®
NATVNEG-6C	Vaginal Negative Control (6 x 0.9 mL)	IVD	Cepheid® GeneXpert®
NATVPOS-6C	Vaginal Positive Control (6 x 0.9 mL)	IVD	Cepheid® GeneXpert®

Healthcare-Associated Infections

Part Number	Product Description	Use*	Platform(s)**
NATCSO-6MC-IVD	C. difficile Negative Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATCDI-6MC-IVD	C. difficile Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATMRSA-6MC-IVD	MRSA Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATMSSE-6MC-IVD	MRSA/SA Negative Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®
NATMSSA-6MC-IVD	S. aureus (MSSA) Positive Control (6 x 0.5 mL)	CE IVD 	Cepheid® GeneXpert®

* RUO = Research Use Only; CE IVD = In Vitro Diagnostic with European CE mark, Notified Body code 2797 applied as applicable

UK Conformity Assessed (UKCA) marking is a conformity mark that indicates conformity with the applicable requirements for products sold within Great Britain, Notified Body code 0086 applied to product labeling as applicable



The Americas:
1.800.274.5487
zepto.customerservice@antylia.com
Europe, Middle East, & Africa:
+44 1480 272279
zepto.EMEA@antylia.com

China:
+86 21 5109 9909
zepto.China@antylia.com
India:
+91 22 6139 4444
zepto.India@antylia.com

Australia & New Zealand:
+61738293996
zepto.AustraliaNewZealand@antylia.com
East Asia:
zepto.EastAsia@antylia.com